

Hidden Isochema spares and accessories: IMI

Part Code	Description	Quantity (per set)
IMI001	Ni Foil Burst Disc	1
IMI002	Proportional Relief Valve (PSV1)	1
IMI003	Replacement PSV1 Spring	1
IMI004	Pneumatically-actuated Diaphragm Valve (Normally Open)	1
IMI005	Pneumatically-actuated Diaphragm Valve (Normally Closed)	1
IMI006	Solenoid Valve	1
IMI007	IMI Reactor Thermocouple Gland	1
IMI008	VCR-8 Male Nut (for reactor)	1
IMI009	Type K Thermocouple (1/16" diameter)	1
IMI010	Valve Diaphragm Kit	1
CAB001	Hidden Isochema RS232 Communication Cable 5m	1
CAB004	RS232 Communication Cable (for Grant waterbaths)	1
GSKT003	VCR-8 Stainless Steel Gasket (e.g. for IMI reactor)	10, 25 or 100
GSKT004	VCR-4 Stainless Steel Gasket	10, 25 or 100
GSKT005	VCR-4 Stainless Steel 2µm Filter Gasket	1