

IGA Accessories

Part Code	Description	Quantity (Per set)
GRAV001	Standard Stainless Steel reactor	1
GRAV002	Quartz Reactor tube (with KF40 fitting)	1
GRAV003	DN-35-CF to KF40 adaptor (with ¼" port)	1
GRAV004	Type K Thermocouple (as External or Sample sensor)	1
GRAV005	PRT (Platinum Resistance Thermometer) Sensor (as External sensor)	1
GRAV006	PRT (Platinum Resistance Thermometer) Sensor (as Sample sensor) - gold plated	1
GRAV007	Stainless Steel Swagelok™ 4BK Series Valve	1
GRAV008	Stem Seal Kit for 4BK Series Valve	1
GRAV009	Gasket Kit for 4BK Series Valve	1
GRAV010	Switched Mode PSU (Power Supply Unit) +5 and ± 15 V	1
GRAV011	Switched Mode PSU (Power Supply Unit) +24 V	1
GRAV012	Hangdown Set (includes: 1 sample chain, 1 tare chain, 4 balance hooks, 4 sample hooks)	1
GRAV013	IGA Tungsten Hangdown (diameter 0.125mm, length 300mm)	1
GRAV014	Quartz Hook (diameter 0.7mm, length 120mm)	1
GRAV015	Thermal Sponge	1