

Hidden Isochema Spares – Sample Holders

Part Code	Description	Quantity (per set)
SB001	Standard Stainless Steel Mesh Sample Holder (42µm aperture mesh, approx. diameter 15mm)	5 (can mix and match with Pyrex Sample Holders)
SB002	Standard Stainless Steel Mesh Sample Holder (20µm aperture mesh, approx. diameter 15mm)	5 (can mix and match with Pyrex Sample Holders)
SB003	Standard Stainless Steel Mesh Flat Base Sample Holder (42µm aperture mesh, approx. height 15mm, approx. diameter 15mm)	1
SB004	Large Stainless Steel Mesh Cylindrical Sample Holder (42µm aperture mesh, approx. height 23mm, approx. diameter 15mm)	1
SB005	Stainless Steel Foil Cylindrical Sample Holder (approx. height 11mm, approx. diameter 11mm)	1
SB006	Large capacity Stainless Steel Mesh Cylindrical Sample Holder (42µm aperture mesh, approx. height 23mm, approx. diameter 20mm)	1
SB007	Quartz Sample Holder with Basal Hole (approx. height 11-12.5mm, approx. diameter 11-12.5mm, weight ≤ 0.3g)	1
SB008	Quartz Sample Holder without Basal Hole (approx. height 11-12.5mm, approx. diameter 11-12.5mm, weight ≤ 0.3g)	1
SB009	Large Quartz Sample Holder (approx. height 20-25mm, approx. diameter 20mm, weight ≤ 1g)	1
SB010	Pyrex Sample Holder (approx. bulb diameter 12mm, approx. hole diameter 5mm, weight 250-300mg)	5 (can mix and match with Standard SS Mesh Sample Holders)
SB011	Pyrex Sample Holder (approx. bulb diameter 12mm, approx. hole diameter 5mm, weight 150-200mg)	5 (can mix and match with Standard SS Mesh Sample Holders)