








## Hidden Isochema Spares – Sample Holders

Part Code	Description	Image
SB001	Standard Stainless Steel Mesh Sample Holder (42µm aperture mesh, approx. diameter 15mm)	A stainless steel mesh sample holder with a conical shape and a hanging wire, shown against a grey background.
SB002	Standard Stainless Steel Mesh Sample Holder (20µm aperture mesh, approx. diameter 15mm)	A stainless steel mesh sample holder with a conical shape and a hanging wire, shown against a grey background.
SB003	Standard Stainless Steel Mesh Flat Base Sample Holder (42µm aperture mesh, approx. height 15mm, approx. diameter 15mm)	A stainless steel mesh sample holder with a flat base and a hanging wire, shown against a blue background.
SB004	Large Stainless Steel Mesh Cylindrical Sample Holder (42µm aperture mesh, approx. height 23mm, approx. diameter 15mm)	A large stainless steel mesh cylindrical sample holder with a hanging wire, shown against a blue background.
SB005	Stainless Steel Foil Cylindrical Sample Holder (approx. height 11mm, approx. diameter 11mm)	A stainless steel foil cylindrical sample holder with a hanging wire, shown against a blue background.
SB006	Large capacity Stainless Steel Mesh Cylindrical Sample Holder (42µm aperture mesh, approx. height 23mm, approx. diameter 20mm)	A large capacity stainless steel mesh cylindrical sample holder with a hanging wire, shown against a blue background.
SB007	Quartz Sample Holder with Basal Hole (approx. height 11-12.5mm, approx. diameter 11-12.5mm, weight ≤ 0.3g)	A clear quartz sample holder with a basal hole and a hanging wire, shown against a grey background.
SB008	Quartz Sample Holder without Basal Hole (approx. height 11-12.5mm, approx. diameter 11-12.5mm, weight ≤ 0.3g)	A clear quartz sample holder without a basal hole and a hanging wire, shown against a red background.

## Hidden Isochema Spares – Sample Holders

SB009	Large Quartz Sample Holder (approx. height 20-25mm, approx. diameter 20mm, weight $\leq 1g$ )	
SB010	Pyrex Sample Holder (approx. bulb diameter 12mm, approx. hole diameter 5mm, weight 250-300mg)	
SB011	Pyrex Sample Holder (approx. bulb diameter 12mm, approx. hole diameter 5mm, weight 150-200mg)	
SB012	XEMIS Large Capacity Bucket (approx. diameter 9mm, height 41mm, weight $\leq 0.2g$ )	
SB013	XEMIS Sample Holder SS Foil (approx. diameter 9mm, height 12.5mm, weight $\leq 0.45g$ )	
SB014	XEMIS Sample Holder SS Mesh (approx. diameter 9mm, height 19mm, weight $\leq 0.040g$ )	
SB015	XEMIS Sample Bucket Loading Tool (approx. diameter 9mm, height 12.5mm)	
SB016	IGA Metallic Blue Anodized Sample Holder (approx. base diameter 21mm, inner diameter 11mm, height 10mm)	

## Hidden Isochema Spares – Sample Holders

SB017	Membrane Sample Bucket Base (includes 6 screws and 6 washers)	 A photograph showing a circular, metallic Membrane Sample Bucket Base with six mounting holes around its perimeter. To the left of the base are six screws and six washers, arranged in a small cluster.
SB018	Membrane Sample Bucket Lid	 A photograph of a circular, metallic Membrane Sample Bucket Lid, which is a shallow, wide-rimmed ring with six mounting holes around its outer edge.