

Hidden Isochema Spares and Accessories: IGA

GRAV001	Standard Stainless Steel reactor	1
GRAV002	Quartz Reactor tube (with KF40 fitting)	1
GRAV003	DN-35-CF to KF40 adaptor (with ¼" port)	1
GRAV004	Type K Thermocouple (as External or Sample sensor)	1
GRAV005	PRT (Platinum Resistance Thermometer) Sensor (as External sensor)	1
GRAV006	PRT (Platinum Resistance Thermometer) Sensor (as Sample sensor) - gold plated	1
GRAV007	Stainless Steel Swagelok™ 4BK Series Valve	1
GRAV008	Stem Seal Kit for 4BK Series Valve	1
GRAV009	Gasket Kit for 4BK Series Valve	1
GRAV010	Switched Mode PSU (Power Supply Unit) +5 and ± 15 V	1
GRAV011	Switched Mode PSU (Power Supply Unit) +24 V	1
GRAV012	Hangdown Set (includes: 1 sample chain, 1 tare chain, 4 balance hooks, 4 sample hooks)	1
GRAV013	IGA Tungsten Hangdown (diameter 0.125mm, length 300mm)	1
GRAV014	Quartz Hook (diameter	1

	0.7mm, length 120mm)	
GRAV015	Thermal Sponge	1
CAB001	Hidden Isochema RS232 Communication Cable 5m	1
CAB002	Hidden Isochema RS232 Communication Cable 3m	1
CAB003	IGA Direct Motor Cable	1
CAB004	RS232 Communication Cable (for Grant waterbaths)	1
GSKT001	DN-35-CF Copper Gasket (for IGA reactor flange)	5 or 100
GSKT002	DN-100-CF Copper Gasket (for IGA balance flange)	3
GSKT004	VCR-4 Stainless Steel Gasket	10, 25 or 100
BOLT001	Standard Reactor Bolt Set (includes M6 × 35 bolts, washers and nuts)	6
BOLT002	Quartz Reactor Bolt Set (includes M6 × 25 bolts and nuts)	6
BOLT003	Dry Loader Sealing Screws	6